



Operating Light(DX)

KDLED700(DX)



LED Shadowless Lamp

1. This is the perfect new mold for the trendy butterfly type.
2. It adopts import German quality arm to rotate and fix well.
3. It has a shadowless compensation for anywhere. That makes the illuminance very well.
4. We use OSRAM lamp 73-piece long life lamp system.
5. Advanced control system for automatic focus.
6. It has three types for only 1 button: Mirror type / Glossy type / normal type.




Detail Pareameter:

Parameter	4 4	5 5	6 6	7 7
Illuminance	100,000	120,000	140,000	160,000
Color temperature	3600-5000K setting			
Focused beam depth	≥900mm			
Chip for life	≥80000 hours			
Working distance	500-2000mm			
patient's head temperature	No increase			
Operator head temperature	< 1 °C			
Point diameter	100-300mm adjustment			
Color rendering index	85≤Ra≤100			
Chip power	200W12V			
Supply voltage	200V ± 10% 50HZ			
Install height	2500-3150mm			

KDLED700(DX)



Details	
Illuminance (LUX)	40,000-190,000
Color Temp(K)	3500-5500
Diameter of spot	20-26cm
Bulb type	OSRAM >60000 hours
Color rendering index(ra)	>98
Lightbeam depeth	160cm/63 inch
Light adjust range	10%-100%
Temp rise(operator head)	<1°C
Input power	AC100-240V, 50/60HZ
Best install height	2.7-3.1m



Three types by 1 button control:

- 1,Endoscope type suitable endoscope operation
- 2,Bright type Illuminance and color temp to be max
- 3,Normal type





한국엠바이오(주) | Tel : 02-3672-7784 / Fax : 02-6951-2110

본 점 : 서울시 강남구 선릉로 577, POSTOWER | 지 점 : 서울시 성동구 성수일로 10 서울숲ITCT지식산업센터

사업부 : 서울시 강동구 천호대로 1099, 정산타워 | 공장 : 경기도 화성시 향남읍 상신리 900-3

© Korea Medical Bio Corporation Inc. ALL RIGHT RESERVED.